

TURBINE

SPORTS BREATHING TECHNOLOGY

SIZE MATTERS.

Your nose is as unique as your finger print so we have created a sizing guide to help you select the right Turbine for you.

Simply print out this page and carefully cut around the dotted line of a sizing guide and follow the steps below to size your nose.

TIP: Generally, as men are either medium or large and women are either small or medium, we recommend that you cut out the medium guide first.

When printing, select '**Actual Size**' under page sizing options to ensure guides are printed to scale.



1 If you're unsure which size to start with, cut out the medium sizing guide.

Then with the coloured side facing away from you, position the sizing tool under your nose.



2 Once in position, tilt the sizing tool slightly upward so you can either feel or see where the black tip is located.

For a good fit, the black tip should be located at the crease line. In this example, **Turbine Medium** would be the right size.



3 If the black tip is above the crease line, you are most likely **Turbine Large**.

If below this crease line, you are most likely **Turbine Small**.

To confirm, cut-out the other size guide and follow step 1 & 2 again.

Once you know your size, buy your Turbine Pack containing 3 devices and follow the fitting instructions inside the box.

Please note that this is a guide only. The external size of your nose is not always representative of the internal anatomy. A pronounced deviation or internal scar tissue caused by an injury may mean that you cannot fully insert the size indicated by the guide. You may need to choose one size smaller than indicated.

AFTER PRINTING, THE SCALE IS CORRECT IF THIS LINE  MEASURES 37MM (1.46")

